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
INSTRUCTIONS FOR USE

Description:
3Delta Gingiva Mask is a light-curing resin for the additive production of gingiva masks for dental models. The material is optimized for wavelengths of 385 and 405 nm in layer thicknesses of 100 µm and can be used with 3D printers of following manufacturers:

- Rapid Shape [D-Series]
- Asiga [PRO4K UV / Max UV]
- Ackuretta [SOL 130]

Composition:
Preparation made of acrylates and initiators.

Application:
Additive Manufacturing of gingiva masks.

Use:
BUILD PROCESS

➤ **Prepare data set:** The orientation and type of support structures can be found in the system supplier's instructions for use.

➤ **Establish variable process parameters:** The appropriate process parameters are selected according to the specifications of the printer system supplier. Ensure that the correct material-specific set of parameters is selected (product, layer thickness, exposure time).

➤ **Mixing materials:** Stir material to the bottom of the bottle. First filling: Pour material into the material reservoir of the unit. Refill: Refill the material as needed and mix in the reservoir with a silicone spatula. Allow the reservoir with material to rest in the unit for 10 min.

➤ **Start build process.**

CLEANING AND DRYING

➤ **Clean part:** Separate the parts (with support structures) from the build platform with a spatula. The component is cleaned under flowing isopropanol for max. 5 min. If necessary, carefully clean the objects with a brush and a spray bottle filled with isopropanol. ➤ **Dry part:** Blow component with compressed air and leave to dry for 30 min at room temperature or in a drying cabinet adjusted to a temperature of approx. 37 °C.

POST-CURING

➤ For post-curing, the components have to be cleaned and cooled to room temperature.

Used printer	Curing Device	Post-curing
Rapidshape	Otoflash	2000 flashes each for top and bottom ¹
Asiga Ackuretta	RS Cure	Parameter set <i>3Delta Gingiva Mask</i> ²
Ackuretta	CURIE	T(10) P(14) D (10) B (on) ³

¹ Curing must be performed under inert gas flooding or oxygen exclusion in order to prevent the formation of the inhibition layer. ² exposure time: 5 min / vacuum 50 mbar / without nitrogen / power 75 % / wavelength 375+415 nm ³ (T) Time / (P) Power / (D) Duty Cycle / (B) Bottom Light

The device used must be maintained and calibrated according to the manufacturer's specifications to ensure mechanical properties and color.

➤ Remove support structure if necessary.

➤ Blow component with compressed air and leave to dry.
➤ POST-PROCESSING COMPLETED

Further processing instructions:
Parts are subsequently processed and finalized using standard dental techniques.

Warnings:
May cause an allergic skin reaction. Causes skin irritation. Avoid contact with unpolymerized material. Wear protective gloves, protective clothing, eye protection. If on skin: Wash with plenty of water and soap. If skin irritation or rash occurs: Get medical advice/attention. Avoid breathing vapours/ spray. Switch on extraction system. The material is harmful to the environment. Avoid release to the environment. Collect spillage. Please always observe the respective valid safety data sheets.Always observe the relevant country-specific health and safety regulations and legal requirements when handling hazardous substances.

Storage instructions / shelf life:
Protect from exposure to light. May polymerize spontaneously. Store between 4 °C and 25°C. Close after use. Do not use after the expiration date. See imprint/label for expiration date.

Disposal instructions:
Dispose of contents/container in accordance with local regulations.

The product was developed for use in dentistry and must be used according to the instructions for use.

	Batch number		Keep away from light
	Item number		Temperature limit
	Manufacturer		Observe directions for use
	Use-by date		Stir before use
	Exclamation mark		Envrionmental Hazard

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