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## INSTRUCTIONS FOR USE en

**Description:**  
3Delta Model Blue / Dark Beige / Sand / Grigio are light-curing resins for the additive production of dental models. The materials are optimized for wavelengths of 385 and 405 nm in layer thicknesses of 50 µm and/or 100 µm and can be used with 3D printers of following manufacturers:

- Rapid Shape [D-Series]
- Asiga [PRO4K UV / Max UV]
- Ackuretta [SOL 130]

**Composition:**  
Preparation made of acrylates and initiators.

**Application:**  
Additive manufacturing of dental models, in particular for prosthetic master models and orthodontic models e.g. for aligner manufacturing or situation models

**Use:**  
**BUILD PROCESS**  
 ↳ **Prepare data set:** The orientation and type of support structures can be found in the system supplier's instructions for use.  
 ↳ **Establish variable process parameters:** The appropriate process parameters are selected according to the specifications of the printer system supplier. Ensure that the correct material-specific set of parameters is selected (product, layer thickness, exposure time).  
 ↳ **Mixing materials:** Stir material to the bottom of the bottle. First filling: Pour material into the material reservoir of the unit. Refill: Refill the material as needed and mix in the reservoir with a silicone spatula. Allow the reservoir with material to rest in the unit for 10 min.  
 ↳ **Start build process.**

**CLEANING AND DRYING**  
 ↳ **Option 1:** Separate part (with support structures) from the build platform with a spatula. Place the part in the RS Wash (Rapid Shape) and clean with isopropanol, selecting the appropriate parameter set. Replace the isopropanol after 20 cycles. Other relevant: → *IfU RS Wash*.  
 ↳ **Option 2:** Separate the parts (with support structures) from the build platform with a spatula. The component is cleaned under flowing isopropanol for max. 5 min. If necessary, carefully clean the objects with a brush and a spray bottle filled with isopropanol.  
 ↳ **Dry part:** Blow component with compressed air and leave to dry for 30 min at room temperature or in a drying cabinet adjusted to a temperature of approx. 37 °C.

**POST-CURING**  
 ↳ For post-curing, the components have to be cleaned and cooled to room temperature.

Used printer	Curing Device	Material	Post-curing
Rapidshape Asiga Ackuretta	Otoflash	3Delta Model -Blue -Dark Beige -Sand -Grigio	2000 flashes each for top and bottom <sup>1</sup>
	RS Cure		Selection of the corresponding material parameter set <sup>2</sup>
Ackuretta	CURIE		(T) 20, (P) 16, (D) 10, (B) on <sup>3</sup>

- <sup>1</sup> Recommendation: Curing under exclusion of oxygen prevents the formation of the inhibition layer.
- <sup>2</sup> exposure time: 10 min / vacuum 50 mbar / without nitrogen / power 100 % / wavelength 375+415 nm
- <sup>3</sup> (T) Time / (P) Power / (D) Duty Cycle / (B) Bottom Light

The device used must be maintained and calibrated according to the manufacturer's specifications to ensure mechanical properties and color.

↳ Remove support structure if necessary.

↳ For post-curing without oxygen connection, clean under flowing conditions with isopropanol for max. 2 min. If necessary, carefully clean the objects with a brush and a spray bottle filled with isopropanol.  
 ↳ Blow component with compressed air and leave to dry.  
**↳ POST-PROCESSING COMPLETED**

**Further processing instructions:**  
Parts are subsequently processed and finalized using standard dental techniques.  
**Warnings:**

May cause an allergic skin reaction. Avoid contact with unpolymerized material. Wear protective gloves, protective clothing, eye protection. If on skin: Wash with plenty of water and soap. If skin irritation or rash occurs: Get medical advice/ attention. Avoid breathing dust/ vapours/ spray. Switch on extraction system. The material is harmful to the environment. Avoid release to the environment. Collect spillage. Please always observe the respective valid safety data sheets. Always observe the relevant country-specific health and safety regulations and legal requirements when handling hazardous substances.

**Storage instructions / shelf life:**  
Protect from exposure to light. May polymerize spontaneously. Store between 4 °C and 25 °C. Close after use. Do not use after the expiration date. See imprint/label for expiration date.

**Disposal instructions:**  
Dispose of contents/container in accordance with local regulations.

The product was developed for use in dentistry and must be used according to the instructions for use.

	Batch number		Keep away from light
	Item number		Temperature limit
	Manufacturer		Observe directions for use
	Use-by date		Stir before use
	Exclamation mark		Environmental Hazard

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