de	Gebrauchsanweisung
en	Instructions for Use
fr	Mode d’emploi
es	Instrucciones de uso
it	Istruzioni per l’uso
pt	Instruções de utilização

3Delta Tray Rose 1000 g

REF 1621620

DeltaMed GmbH

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INSTRUCTIONS FOR USEen

Description:
3Delta Tray Rose is a light-curing resin for the additive production of dental impression trays. The material is optimized for wavelengths of 385 and 405 nm in layer thicknesses of 150 µm and can be used with 3D printers of follow-ing manufacturers:

- Rapid Shape [D-Series, PRO-Series]
- Asiga [PRO4K UV / MAX UV]
- Elegoo [Saturn 4 Ultra]

Composition:
Preparation made of acrylates and phosphinate.

Indication:
Additive manufacturing of dental impression trays.

Contraindications:
Patient and user: Product should not be used in patients or by users who are known to be allergic to any of the ingredients. Otherwise the patient group is not restricted.

Adverse reactions:
Patient and user: Product can cause allergic reactions.
User: Observe warnings.

Interactions:
Interactions are not known.

Use:
BUILD PROCESS
➤ **Prepare data set:** The orientation and type of support structures can be found in the system supplier's instructions for use.
➤ **Establish variable process parameters:** The appropriate process param-eters are selected according to the specifications of the printer system sup-plier. Ensure that the correct material-specific set of parameters is selected (product, layer thickness, exposure time).
➤ **Mixing materials:** Stir material to the bottom of the bottle.
First filling: Pour material into the material reservoir of the unit.
Refill: Refill the material as needed and mix in the reservoir with a silicone spatula. Allow the reservoir with material to rest in the unit for 10 min.
➤ **Start build process.**

CLEANING AND DRYING
➤ **Clean part:** Separate the parts (with support structures) from the build plat-form with a spatula. The component is cleaned under flowing isopropanol for **max. 5 min.** If necessary, carefully clean the objects with a brush and a spray bottle filled with isopropanol. Cleaning with RS Wash / PROWASH (Rapidshape) devices is conducted by selecting the corresponding material parameter set (120s pre-cleaning / 65s cleaning / 90s drying / Nr. Of uses: 10).
➤ **Dry part:** Blow component with compressed air and leave to dry for 30 min at room temperature or in a drying cabinet adjusted to a temperature of approx. 37 °C.

POST-CURING
➤ For post-curing, the components have to be cleaned and cooled to room temperature.

Used printer	Curing Device	Post-curing
Rapidshape Asiga elegoo	Otoflash	2000 flashes each for top and bottom / Σ 4000 flashes ¹
	RS Cure / PROCURE	Parameter set <i>3Delta Tray Rose</i> ²

¹ Recommendation: Curing under exclusion of oxygen prevents the formation of the inhibition layer.
² exposure time: 6 min / vacuum 50 mbar / without nitrogen / power 100 % / wavelength 375+415 nm

Ensure that the unit used is serviced and calibrated in accordance with the manufacturer's instructions in order to preserve its mechanical properties, color and biocompatibility.

➤ Remove support structure if necessary.

➤ For post-curing without oxygen connection, clean under flowing condi-tions with isopropanol for max. 2 min. If necessary, carefully clean the ob-jects with a brush and a spray bottle filled with isopropanol.
➤ Blow component with compressed air and leave to dry.
➤ **POST-PROCESSING COMPLETED**

Further processing instructions:
Parts are subsequently processed and finalized using standard dental tech-niques.

Warnings:
Any deviation from the processing process described here can jeopardize the biocompatibility, functionality and safety of the final product. Contains skin sensitizing substances (acrylates, phosphinate). Please always observe the re-spective valid safety data sheets. Applies to the uncured material: May cause

an allergic skin reaction. Avoid contact with unpolymerized material. Wear protective gloves / protective clothing / eye protection. If on skin: Wash with plenty of soap and water. Avoid breathing dust / vapours / spray. Switch on extraction system. If skin irritation or rash occurs: Get medical advice / atten-tion. The unpolymerized material is harmful to the environment. Avoid re-lease to the environment. Collect spillage. Always observe the relevant coun-try-specific health and safety regulations and legal requirements when han-dling hazardous substances.

Storage instructions / shelf life:
Protect from exposure to light. May polymerize spontaneously. Store be-tween 4 °C and 25 °C. Close after use. Do not use after the expiration date. See imprint/label for expiration date.

Disposal instructions:
Dispose of contents/container in accordance with local regulations.

Rx Important product information
Federal (U.S.) law restricts the finished device to sale by or on the order of a dentist.

The product was developed for use in dentistry and must be used according to the instructions for use. Serious incidents must be reported to DeltaMed GmbH and the responsible authority.

	Batch number		Keep away from light
	Item number		Temperature limit
	Manufacturer		Consult instructions for use
	Use-by date		Stir before use
	Exclamation mark		Environmental Hazard
	Medical Device		

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